

Tachy Episode (Archive)

Tachy Episode 3 of 3	Date/Time: 29 Jun 2008 09:50	Type: VT (Monitor)	Episode Duration: 0:18
-----------------------------	-------------------------------------	---------------------------	-------------------------------

Therapy: (Monitor Only)	Results: Below Rate Detection (CL 610 ms)	Alerts: 1 SVT discriminators disagree
---------------------------------------	---	---

	<p style="text-align: center;">Diagnosis: VT (Monitor)</p> <p style="text-align: center;">Time to Diagnosis: 8.00 sec CL: 320 ms/187 bpm</p> <p>VT Diagnosis Criteria: 2 of 3</p> <p>Morphology:</p> <p style="padding-left: 20px;">Programmed: On, $\geq 60\%$ is a match, ≥ 5 of 8 matches indicate SVT Template Updated: 22 Jun 2008 06:54</p> <p style="padding-left: 20px;">Measured:</p> <p style="padding-left: 40px;">Min Match Score: N/A</p> <p style="padding-left: 40px;">Max Non-Match Score: N/A</p> <p style="padding-left: 40px;">No. Template Matches: N/A (VT Indicated)</p> <p>Interval Stability:</p> <p style="padding-left: 20px;">Programmed: On, ≥ 60 ms indicates SVT (SIH Passive, ≥ 2 intervals indicates SVT)</p> <p style="padding-left: 20px;">Measured:</p> <p style="padding-left: 40px;">Stability Delta: 95 ms (SVT Indicated)</p> <p>Sudden Onset:</p> <p style="padding-left: 20px;">Programmed: On, $<20\%$ indicates SVT</p> <p style="padding-left: 20px;">Measured Max Delta: 30 % (VT Indicated)</p>
--	---

SVT Criteria Statistics

Number of SVT Diagnoses during this Episode: 0

Measured SVT Discriminator values for Episode:

Min / Max Stability Delta: 95 ms / 95 ms

Max Onset Delta: 30 %



Tachy Episode (Archive)

Tachy Episode 3 of 3		Type: VT (Monitor)		Episode Duration: 0:18	
Therapy:		Results:		Alerts: 1	
(Monitor Only)		Below Rate Detection (CL 610 ms)		SVT discriminators disagree	

<p style="text-align: center;">Diagnosis: VT (Monitor)</p> <p style="text-align: center;">Time to Diagnosis: 8.00 sec CL: 320 ms/187 bpm</p> <p>VT Diagnosis Criteria: 2 of 3</p> <p>Morphology:</p> <p>Programmed: On, ≥ 60% is a match, ≥ 5 of 8 matches indicate SVT Template Updated: 22 Jun 2008 06:54</p> <p>Measured:</p> <p>Min Match Score: N/A</p> <p>Max Non-Match Score: N/A</p> <p>No. Template Matches: N/A (VT Indicated)</p> <p>Interval Stability:</p> <p>Programmed: On, ≥ 60 ms indicates SVT (SIH Passive, ≥ 2 intervals indicates SVT)</p> <p>Measured:</p> <p>Stability Delta: 95 ms (SVT Indicated)</p> <p>Sudden Onset:</p> <p>Programmed: On, <20 % indicates SVT</p> <p>Measured Max Delta: 30 % (VT Indicated)</p>	
<p style="text-align: center;">SVT Criteria Statistics</p> <p>Number of SVT Diagnoses during this Episode: 0</p> <p>Measured SVT Discriminator values for Episode:</p> <p>Min / Max Stability Delta: 95 ms / 95 ms</p> <p>Max Onset Delta: 30 %</p>	

