



Tachy Episode (Archive)

Tachy Episode 2 of 2		Type: VF (Therapy Delivered)	Episode Duration: 0:24
Diagnosis: VF (SVT Discrimination is Off in this zone) Time to Diagnosis: 2.50 sec CL: 165 ms/363 bpm		Alerts: 1 At least one shock unsuccessful	
Therapy: Defib 15 J / 537V Defib 36 J / 830V		Results: VF Below Rate Detection (CL 355 ms)	
		High Voltage Therapy Details First Charge Time: 2.7 s Last Charge Time: 9.4 s Last Lead Impedance: 38 Ω Delivered PW: +4.9 ms, -4.9 ms	

EGM 1: A Sense / Pace (2.2 mm/mV, ± 3.0 mV) EGM 2: V Sense / Pace (1.5 mm/mV, ± 8.9 mV) Sweep Speed: 25.0 mm/s

Atlas™+ DR V-243

Tachy Episode (Archive)

Atlas™+ DR V-243
Merlin™ PCS

ST. JUDE MEDICAL™

Tachy Episode page 2 of 2